

# Hot AXIOS™


Stent and Electrocautery  
Enhanced Delivery System



*The first stent CE Marked for endoscopic  
ultrasound guided transluminal therapy*

Cautery-enabled access and delivery catheter  
with the preloaded therapeutic Hot AXIOS Stent  
for an exchange-free procedure.


- Proprietary one-step combined diathermic ring and cut-wire provides easy access into target tissue
- Perpendicular flanges secure tissue layers and help to prevent migration
- Stent creates anastomotic conduit between two lumens
- Large diameter fully covered lumen apposing stent enables drainage and passage of the endoscope for additional therapeutic procedures



# Hot AXIOS™ Stent and Electrocautery Enhanced Delivery System

Indicated for use to facilitate transgastric or transduodenal endoscopic drainage of:

- A pancreatic pseudocyst or a walled-off necrosis with  $\geq 70\%$  fluid content
- Gallbladder in patients with acute cholecystitis who are at high risk or unsuitable for surgery
- The bile duct after failed ERCP in patients with biliary obstruction due to a malignant stricture




## Ordering Information

### Hot AXIOS™

Stent and Electrocautery Enhanced Delivery System

Order Number	Lumen Diameter (mm)	Saddle Length (mm)	Flange Diameter (mm)	Catheter Outer Diameter (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)
M00553520	6	8	14	9	138	144
M00553530	8	8	17	9	138	144
M00553540	10	10	21	10.8	138	144
M00553550	15	10	24	10.8	138	144
M00553570*	15	15	24	10.8	138	144
M00553560	20	10	29	10.8	138	144



Compatible with Therapeutic EUS scopes with working channel 3.7mm or more.

\* The 15 x 15 HOT AXIOS Stent and Electrocautery-Enhanced Delivery System is indicated for use to facilitate transgastric or transduodenal endoscopic drainage of a pancreatic pseudocyst or a walled-off necrosis with  $\geq 70\%$  fluid content.

**WARNING:** The safety and effectiveness of the 20mm Stent for drainage of the Gallbladder and Bile Duct has not been established.

Contact your local Boston Scientific representative for more information.

**Boston Scientific**  
Advancing science for life™

All trademarks are the property of their respective owners.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France. All photographs taken by Boston Scientific.

ENDO-458506-AC. Produced by Gosling.

www.bostonscientific.eu

2021 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved. DINEND2374EC