


BD Instaflash™ Technology - Designed to Increase First-stick Success



- A successful first-stick is important for both the user and the patient, especially when accessing small and compromised veins
- BD Instaflash™ Needle Technology is designed to reduce hit-and-miss insertion by immediately confirming vessel entry through a unique notch in the needle shaft

Products Manufactured with BD Instaflash™

- BD Insyte™
- BD Insyte-N™
- BD Insyte™ Autoguard™
- BD Insyte-N™ Autoguard™
- BD Insyte™ Autoguard™ BC
- BD Nexiva™
- BD Nexiva Diffusics™
- BD Insyte-A™ Arterial Catheters
- BD Venflon™ Pro Safety
- BD Cathena™
- BD Neoflon™ Pro Safety