



REAGENT STRIPS FOR URINALYSIS

Production Profile

DIRUI

REAGENT STRIPS FOR URINALYSIS

Production Profile

Urinalysis Strip Description	Test Parameters														Packing
	Leucocytes	Nitrite	Urobilinogen	Protein	pH	Blood	SG	Ketone	Billirubin	Glucose	Ascorbic acid	Microalbumin	Creatinine	Calcium	
FUS-10	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙					
FUS-11	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙				
FUS-11MA	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙		⊙			
FUS-12MA	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙			
FUS-13Cr	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙		
FUS-14Ca	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	
H10-800	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙					
H11-800	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙				
H11-800MA	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙		⊙			
H12-800MA	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙			
H13-Cr	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙		
H14-Ca	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	⊙	

100 strips/bottle

Technical Advantage:

- Having the ability of resisting the interference.
- Suitable for full-automatic analyzers and raising the efficiency
- Avoiding the phenomenon of misdiagnosis and improving the reliability of the test result.
- High sensitivity, which can test the minium content of pathology and help to patient's early diagnosis and treatment.
- Microalbumin-creatinine ratio helps laboratories reflect early renal damage such as diabetic nephropathy
- Specialization, not to react with other non-reactive substance and make the test result more accurate



Certified to
ISO 9001:2008 and ISO 13485:2003

DIRUI INDUSTRIAL CO., LTD.

3333 Yiju Street, New & High Tech. Development Zone
Changchun, Jilin 130103, P.R. China

Tel: +86(431)81935329 85100409

Fax: +86(431)85172581 85083741

E-mail: dirui@dirui.com.cn Http://www.dirui.com.cn

• Specifications subject to change without notice.

20170413

DIRUI