



# BD Totalys™ SlidePrep

Efficiency and confidence in slide processing

- Processes BD SurePath™ Liquid-based Pap test and Non-Gynaecological samples
- Integrates slide creation and staining from an enriched cell pellet
- Allows for random sample loading
- Fully integrates with the BD Totalys™ MultiProcessor to maintain chain-of-custody




# Help improve your lab's productivity

## A new level of confidence through automation

- Consistent results with staining in individual chambers using fresh reagents for every sample
- Accurate reagent pipetting and precise timing for each sample reduces process variability
- Automatic barcode scanning ensures sample match and chain of custody

## Efficient throughput

- Complete a single run of 48 slides in approximately 1 hour<sup>2</sup>
- High throughput capability: up to 336 stained slides produced from cell pellets in an 8-hour shift<sup>2</sup>
- Minimal hands-on time allows staff to perform other tasks



BD Totalys™ SlidePrep

## Set the precedent

Contact your BD representative today.

Specimen handling	
Specimen type	Enriched cell pellets and non-gynaecological cell pellets
On-board sample capacity	Up to 48 samples
Positive sample identification	Barcodes on all centrifugation tubes and cytology slides ensure chain of custody throughout processing

Dimensions	
Height	86.4 cm (34 in)
Width	182.9 cm (72 in) (includes bulk reagent station, disposable tip waste container and clearance)
Depth	62.2 cm (24.5 in)
Weight	124.7 kg (275 lbs)
Separate bulk reagent station	35.6 cm H x 55.9 cm W x 35.6 cm D (14 in H x 22 in W x 14 in D)
Separate waste station dimensions	40.6 cm H x 33 cm W x 58.4 cm D (16 in H x 13 in W x 23 in D)

Environmental requirements	
Temperature	18–30 °C (64–86 °F)
Humidity	20–80% RH, non condensing
Altitude	Up to 2,000 m (6,562 ft)

Electrical requirements	
Input voltage	100–240 VAC
Rated current	12 A
Input line frequency	50–60 Hz
Power	1,000 Watts

1. BD Totalys™ SlidePrep Instrument User's Manual [500005570]. Franklin Lakes, NJ: Becton, Dickinson and Company.  
2. Internal BD Document, PD-08129-67.

BD - Europe, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland

[bd.com](http://bd.com)



BD SurePath™, BD Totalys™ Multiprocessor, BD Totalys™ SlidePrep are bearing a CE mark. BD, the BD logo, BD SurePath™, BD Totalys™, are trademarks of Becton, Dickinson and Company. All other trademarks are from their respective owners. © 2022 BD. All rights reserved. BD-59833

