


## BD Vialon™ Biomaterial

- BD Vialon™ Biomaterial is a proprietary polyurethane, developed especially for vascular access.

BD Vialon™ is clinically proven to

- Increase the ease of penetration<sup>1</sup>
- Reduce the risk of infiltration<sup>2</sup>
- Reduce the risk of phlebitis thereby reducing costs and increasing indwell times<sup>1,3,4,5</sup>
- Exhibit greater kink resistance ensuring consistent flow rates<sup>1,3</sup>



## Products Manufactured with BD Vialon™

- |                            |                        |                          |
|----------------------------|------------------------|--------------------------|
| • BD Insyte™               | • BD Insyte-W™         | • BD Saf-T-Intima™       |
| • BD Insyte-N™ Autoguard™  | • BD Insyte-N™         | • BD Venflon™ Pro        |
| • BD Insyte™ Autoguard™    | • BD Nexiva™           | • BD Venflon™ Pro Safety |
| • BD Insyte™ Autoguard™ BC | • BD Nexiva Diffusics™ |                          |

1. Gaukrager PB, Roberts JG, Manners TA; Infusion thrombophlebitis: A prospective comparison of 645 Vialon™ and Teflon® cannulae in anaesthetic and post-operative use. *Anaesth. Intens. Care* 1988; 16:265 - 271
2. Stanley MD, Meister E, Fuschuber K; Infiltration during intravenous therapy in neonates: Comparison of Teflon® and Vialon™ catheters. *South Med J* 1992; 85(9):883 - 886.
3. Jacquot C, Fauvage B, Bru JP, Croize J, Calop J; Effect of type of material on thrombophlebitis risk with peripheral venous catheters. *Ann Fr Anesth Réanim* 1989; 8:3 - 7.
4. Maki DG, Ringer M; Risk factors for infusion-related phlebitis with small peripheral venous catheters: a randomised controlled trial. *Ann. Intren. Med* 1991; 114:845 - 854.
5. McKee JM, Shell JA, Warren TA, Campbell VP; Complications of intravenous therapy; a randomised prospective study - Vialon™ vs Teflon®. *J Intrav. Nurs*, 1989; 12(5): 288 - 295